NERI®

NOTA INFORMATIVA cod. 392036 - mod. Flexto Blu Vinile

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NERI®

INFORMATIVE NOTE cod. 392036 - mod. Flexto Blu Vinile

Read this briefing note through carefully prior to use and prior to any maintenance work. The purpose of the information contained herein is to assist and direct the user in choosing and using the equipment. Neither the manufacturer nor the distributor can be held liable for an incorrect use of the equipment. This briefing note must be kept for the duration of the equipment. Glove model: **cod. 392036 - mod. Flexto Blu Vinile.**

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CE2777

The CE 2777 mark means that these gloves are personal protective equipment conforming to the essential health and safety requirements set forth n the (EU) Regulation 2016/425 and have been certified by the notified body: **SATRA Technology Europe Ltd, Bracetown Business Park, Clonee, Dublin, D15YN2P, Ireland n° 2777.** In addition, the equipment is subject to the control of CE quality guarantee of the finish product by the notified body: **SATRA Technology Europe Ltd, Bracetown Business Park, Clonee, Dublin, D15YN2P, Ireland n° 2777.**

SIZES

OUR SIZE REFERENCES	SIZES CONTEMPLATED EN 420:2003+A1:2009
S	6
M	7
L	8
XL	9

The dexterity level of the fingers is equivalent to 5 out of 5.

The gloves meet the safety characteristics only if the size is right and worn properly.

RANGE OF USE

The glove is designed to protect the user against contact with chemical products and micro-organisms (intended as fungus, bacteria and viruses).

PICTOGRAMS AND PROTECTION LEVELS

 EN ISO 374-1:2016 TYPE B	 EN ISO 374-5:2016
KPT	VIRUS

PROTECTION AGAINST CHEMICALS	PROTECTION AGAINST MICRO-ORGANISMS Resistance to air and water penetration using the method referred to in EN 374-2:2014 standards
-------------------------------------	--

CHEMICAL RESISTANCE TABLE according to EN ISO 374-1: 2016

Chemical	Permeation index [0 - 6]	Time of Passage [min]	Degradation index
K Sodium hydroxide 40%	5	240	-1.3%
P Hydrogen peroxide 30%	5	240	0.7%
T Formaldeide 37%	6	480	-11.3%

Time of passage. The time a chemical takes to pass through the glove material. The values may differ when the chemical substance is used in a compound.

Do not use the gloves for handling the above-mentioned chemical products for longer than the times stated above.

 CONFORME A: RÉGLEMENTATION (UE) N. 10/2011 - RÉGLEMENTATION (CE) N. 1935/2004 - DIRECTIVES 85/572/EEC AND SUBSEQUENT AMENDMENTS.			
Ethanol 50%	Acetic acid 3%	Vegetable oil	Polyphenyl oxide (MPPD)
Yes	Yes	No	No
2 h / 40°C	1 h / 40°C	-	-

Since the performance levels are based on lab test results, they may differ from the actual conditions at the workplace. Consequently it is advisable to conduct a preliminary test of the gloves to see if they are suitable for their intended use.

The protective characteristics refer to new equipment, in good condition and that has never been subject to cleaning treatments.

WARNINGS

Prior to use, make sure the gloves are in good condition by which we mean no cuts, holes or cracks, etc. If the gloves are not in perfect condition, replace them immediately. The gloves must be used solely for the risks contemplated in this briefing note.

Movements, tears, frictions and degradation caused by contact with the chemical substance may significantly reduce the actual time of use. In the case of corrosive agents, degradation may be the most important factor to consider when choosing the gloves.

STORING

The gloves are packed in cardboard boxes.

Keep the gloves in their original packaging in a clean, dry place, away from sources of heat and direct sunlight. The user must always inspect the gloves to ensure they are sound prior to using them. The glove must be used, provided the packaging is perfectly intact and kept in the preservation conditions stated on the packaging itself, **within the indicated expiry date and no later.** Since these are disposable gloves to be used just the once, they must be replaced at the end of their use – they cannot be used a second time.

NOTES

Used gloves can be contaminated with chemicals or biological substances: dispose of them in compliance with the local regulations in force on this subject (dump, incinerator). The presence of any allergenic substances is not known to the manufacturer at this time: please report any cases found of hypersensitivity or allergic reactions. In the case of manufacturing defects this equipment will be replaced. *** In the event of discrepancies between the different translations, the Italian version will exclusively apply.** For more information, contact:

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